



QuadraLite™ anaesthetic face masks

A truly unique anaesthetic face mask



Anaesthesia ▪ Anaesthetic Face Masks



QuadraLite™ anaesthetic face masks

A truly unique anaesthetic face mask

Unlike traditional face masks, QuadraLite does not use an inflatable cushion. Innovative moulding techniques have enabled the combination of two ideally suited materials to create a face mask with an anatomical shape and soft seal. This provides a comfortable superior airtight fit with a clear shell to enhance patient face visibility.

Four sizes for all patients

The QuadraLite range has only four sizes instead of the traditional seven, providing a solution for all patients due to its unique anatomical shape and requirement for minimal downward pressure to achieve an airtight seal.



The greener, safer choice

QuadraLite is manufactured from polypropylene and thermoplastic elastomer, contains no PVC and is phthalate-free.

Transparent shell

For excellent visibility

Anatomical shape

Provides a precise fit combined with lower dead space

Soft seal

A comfortable airtight fit with flexibility for different face sizes

Colour-coding

Matches other Intersurgical products for easy selection






www.intersurgical.com/info/quadralite

Video available at www.intersurgical.com





Our quick sizing guide for comparison

Patient	Age	Weight (Kg)	Traditional mask size	QuadraLite™ masks
Neonatal	Preterm	<2	0	 Code: 7191000
	Term	2–4		
Infant	3 months	6	1	
	6 months	8		
	9 months	9		
Paediatric	1 year	10	2	
	2 years	12		
	3 years	14		
	4 years	16		
	5 years	18		
Small adult	6 years	20	3	 Code: 7193000
	7 years	23		
	8 years	26		
	9 years	28		
	10 years	32		
	11 years	35		
	12 years	39		
Medium adult	13 years	44	4	
	14 years	50		
	15 years	54		
	16 years	58		
Large adult	16+ years	60–120	5	 Code: 7194000
Extra large adult	16+ years	>120	6	



Code	Description	Box Qty.
7194000	QuadraLite, anaesthetic face mask, size 4, large adult, green seal, 22F (equivalent to traditional mask sizes 5–6)	35
7193000	QuadraLite, anaesthetic face mask, size 3, medium adult, yellow seal, 22F (equivalent to traditional mask sizes 3–5)	35
7192000	QuadraLite, anaesthetic face mask, size 2, small adult, white seal, 22F (equivalent to traditional mask sizes 2–3)	25
7191000	QuadraLite, anaesthetic face mask, size 1, infant, grey seal, no hook ring, 15M (equivalent to traditional mask sizes 0–1)	30
7089000	Anaesthetic face mask hook ring for both QuadraLite and EcoMask with 22mm port	150

Make an enquiry

Lower environmental impact product



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
 T: +44 (0)118 965 6300 F: +44 (0)118 965 6356 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 9001:2015, ISO 13485:2016 and ISO 14001:2015

Please think before you print
 Save energy and paper.
 If you must print this information sheet please print it double sided.

IS5.24 ▪ Issue 7 07.19