



## Fabian HFO and fabian evolution neonatal ventilator range


Breathing systems suitable for use with acutronic fabian ventilators



**Critical Care** ▪ Breathing Systems


## Fabian HFO™ and fabian evolution™ neonatal ventilator range

A comprehensive range of breathing systems suitable for use with the Acutronic® fabian HFO™ and the fabian evolution™ ventilators, allowing for conventional and high frequency ventilation choices. Available manufactured from Flextube™ or Smoothbore tubing, the range features dual and single heated wire breathing systems, infant nasal CPAP systems as well as a transport option.



High Frequency Oscillation (HFO) adaptor kits are also available allowing the clinician to convert these breathing systems and the transport breathing system to HFO ventilation if required.

We have included an nFlow™ heated wire breathing system and an nFlow adaptor kit to convert any of our featured breathing systems to fit the fabian HFO™ and fabian evolution™ ventilators. With a new detachable monitoring line allowing for out-of-incubator care we now offer a complete care solution for infant nasal CPAP.



Make an enquiry

Code	Description	Box Qty.
4651000	HFO adaptor kit to modify fabian breathing systems to HFO ventilation	20
4650000	10mm Flextube transport breathing system for use with Fabian® ventilators with monitoring line, 15M machine connections and swivel Y-piece. ≥ 1.2m	15
4652000	10mm Flextube transport breathing system for use with Fabian® ventilators with monitoring line, 15M machine connections and swivel Y-piece. ≥ 1.6m	10
4072000	10mm Flextube dual heated wire breathing system for use with Fabian® ventilators with monitor line and extra limb. ≥ 1.6m	12
4071000	10mm Flextube single heated wire breathing system for use with Fabian® ventilators with water trap, monitor line and extra limb. ≥ 1.6m	12
5072000	10mm Smoothbore dual heated wire breathing system for use with Fabian® ventilators with monitor line and extra limb. ≥ 1.2m	20
5071000	10mm Smoothbore single heated wire breathing system for use with Fabian® ventilators with water trap, monitor line and extra limb. ≥ 1.2m	20

Acutronic and Fabian are registered trademarks of Acutronic Medical Systems AG. Fabian HFO and fabian evolution are trademarks of Acutronic Medical Systems AG.

Illustrated

IS1.38 ▪ Issue 5 09.23



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK  
 T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com



The manufacturer Intersurgical Ltd is certified to ISO 14001:2015, ISO 9001:2015, ISO 13485:2016 and MDSAP

**Please think before you print**  
 Save energy and paper.  
 If you must print this information sheet please print it double sided.